

pH Balanced
Compatible with
Amalgam
Separators



BaseVac™ | NeutraClean

A Super Powerful Cleaner

A concentrated cleaner, NeutraClean helps dissolve organic debris, blood, and also sugars found in fluoride gels and prophy paste that can accumulate daily in vacuum tubings and suction lines. Breaking down organic debris and regular cleaning and flushing of the dental lines helps move the sediments such as amalgam and prophy pumice in the suction lines to the amalgam separator for recycling. The cleaning action of the atomizer also helps keep HVE valves from getting stuck (hard to turn) from blood and sugars crystalizing in the valves.

NeutraClean is a super-concentrated, dual action cleaner. It safely and effectively cleans and removes debris in filters, lines and components to enhance suction performance and minimize wear on vacuum pumps. The unique neutral pH balanced formula is also compatible with amalgam separator devices.

Note: no measuring cups or powders to mix, it's simple and easy to use - can save 5 min per operator, since the Atomizer can be left to complete cleaning while dental team can perform other tasks!

Key Features Of NeutraClean

- Dislodges and removes organic (i.e. blood and saliva) and inorganic (i.e. amalgam and fluoride gel) debris.
- Phosphate-free formula.
- Formulated with non-foaming agents to protect pump parts and prolong the life of your vacuum pump.
- Deodorizes with a clean Mountain-Fresh scent.
- Super-concentrated formula - 3/4 ounce makes full atomizer of solution (cleans 3 operatories).

Promotion - Q4 - 2018

Purchase 1 case of 4 bottles of solutions and get free Atomizer (\$100 value) or \$30 off NeutraClean Intro Kit

\$1.50 per day - 3 operator clinic