T-Quick Autoclave —••

The Only Manual Autoclave with Closed Door Drying

Post Covid-19: Chairside sterilization is no longer recommended due to aerosols!

The T-Quick autoclave provides:

- * Faster cycles that deliver dry instruments and no wicking
- * Closed door drying
- * More partial loads when needed vs. full loads



Features:

- T-Quick allows you to easily track each phase of the sterilization process
- Automatic safety shutoff prevents overheating
- Water detector with automatic shutoff
- HEPA air filter provides sterile, bacteria-free air for drying.
- Aligns with AAMI guidelines regarding sterilization cycles.



Sterilization Programs

lcon	Cycle Program	Load Type	Cycle Time (min.)*
℅	Unwrapped 134	Unwrapped instruments	11 Minutes
℅	Wrapped 134	Wrapped instruments	15 Minutes
M	Unwrapped 121	Unwrapped delicate instruments	25 Minutes
	Wrapped 121	Wrapped delicate instruments	25 Minutes
*	Unwrapped	Single Instruments	9 Minutes

* Cycle times are based on:1 kg standard loads (not including drying), Warm chamber (and subject to load size)

Panel Features



1 Cycle Status Indicator

- 2 Pressure Indicator
- **3** Power Switch
- 4 Sterilization Timer
- 5 Water Pump Switch
- 6 Drying and Sterilization Timer
- 7 Night Mode Switch
- 8 Temperature Selector

Technical Specifications:

- Chamber Dimensions W x L (mm) 250x467
- Chamber Volume 23L / 6.1gal
- Standard Cassettes Capacity 3 Full
- No. of Trays (4)
- Tray Dimensions WxHxD (mm) 168x20x414
- Voltage(V) / Freq.(Hz) 230V / (50/60 Hz)
- Current(A) / Power(W) 9.6A / 2200W
- Overall Dimensions WxHxD (mm) 508x362x550







Tuttnauer USA Co. Ltd. 25 Power Drive, Hauppauge, NY 11788 Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720 E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com



© 2020 TuttnauerUSA DLIT0010 Rev A 0120